

Bone Fracture Detection via X-Ray Image Processing

Zhoubin Maneshi

11.04.2025






Overview

- **Dataset:** a collection of radiographic images ($\text{IMG_SIZE}=(224, 224)$, in greyscale) covering all anatomical body regions (limbs, hips, knees, ...)
- **Categories:** train, test, and validation
- **Subcategories:** fractured and non-fractured
- **Challenges:**
 - Dataset contains corrupted images (removed via cleaning)
 - Random file names (unimportant when using **ImageDataGenerator**)
 - Multiple copies of the same image
 - Rotated images

Sample Image Grid



Image Distribution across Splits





Preprocessing and Augmentation

- **Constants:** `IMG_SIZE=(224,224)`, `BATCH_SIZE=32`
- **Random Transformations:** rotation, shift, zoom, flips
- **Normalization:** rescales pixel values to `[0, 1]` for better Neural Network performance
- **ImageDataGenerator** is applied for training, validation, and testing a binary classification task
 - **shuffle=False** for test data ensures prediction
 - Class labels are inferred from subfolder names (frac, non-frac)



Build a Model

- **Architecture:** sequential
- **Convolutional Layers:** 32 filters, ReLU activation, 3x3 Kernel size, 224x224x3 input shape
- **Pooling Layers:** MaxPooling(2, 2) downsamples feature maps
- **Flattenning:** to create 1D dense layers
- **Dropout (0.5)** drops 50% of neurons during training
- **Dense** sets up fully connected layers
- **Final Layer:** 1 neuron (binary classification) and sigmoid activation



Model Summary

Model: "sequential"

Layer (type)	Output Shape	Param #
=====		
conv2d (Conv2D)	(None, 222, 222, 32)	896
max_pooling2d (MaxPooling2D)	(None, 111, 111, 32)	0
conv2d_1 (Conv2D)	(None, 109, 109, 64)	18496
max_pooling2d_1 (MaxPooling2D)	(None, 54, 54, 64)	0
flatten (Flatten)	(None, 186624)	0
dropout (Dropout)	(None, 186624)	0
dense (Dense)	(None, 128)	23888000
dense_1 (Dense)	(None, 1)	129



Training Phase

289/289 [=====] - 631s 2s/step - loss: 0.6680 - accuracy: 0.6179 - val_loss: 0.5920 - val_accuracy: 0.7358
Epoch 2/10

289/289 [=====] - 486s 2s/step - loss: 0.5340 - accuracy: 0.7302 - val_loss: 0.5336 - val_accuracy: 0.7491
Epoch 3/10

289/289 [=====] - 477s 2s/step - loss: 0.4506 - accuracy: 0.7911 - val_loss: 0.3950 - val_accuracy: 0.8359
Epoch 4/10

289/289 [=====] - 424s 1s/step - loss: 0.3823 - accuracy: 0.8274 - val_loss: 0.3719 - val_accuracy: 0.8601
Epoch 5/10

289/289 [=====] - 380s 1s/step - loss: 0.3069 - accuracy: 0.8751 - val_loss: 0.3493 - val_accuracy: 0.8770
Epoch 6/10

289/289 [=====] - 399s 1s/step - loss: 0.2753 - accuracy: 0.8867 - val_loss: 0.2845 - val_accuracy: 0.8854
Epoch 7/10



289/289 [=====] - 401s 1s/step - loss: 0.2343 - accuracy: 0.9090 - val_loss: 0.2563 - val_accuracy: 0.9023
Epoch 8/10

289/289 [=====] - 433s 1s/step - loss: 0.2157 - accuracy: 0.9200 - val_loss: 0.2702 - val_accuracy: 0.9083
Epoch 9/10

289/289 [=====] - 414s 1s/step - loss: 0.1926 - accuracy: 0.9287 - val_loss: 0.2419 - val_accuracy: 0.9156
Epoch 10/10

289/289 [=====] - 398s 1s/step - loss: 0.1764 - accuracy: 0.9332 - val_loss: 0.2196 - val_accuracy: 0.9276

Accuracy and Loss




16/16 [=====] - 4s 248ms/step - loss: 0.2476 - accuracy: 0.9190

Test Accuracy: 91.90%

16/16 [=====] - 4s 214ms/step

	precision	recall	f1-score	support
not_fractured	0.91	0.92	0.91	238
fractured	0.93	0.92	0.92	268
accuracy			0.92	506
macro avg	0.92	0.92	0.92	506
weighted avg	0.92	0.92	0.92	506

Confusion Matrix





Further Improvement

- **Use of Hashing to remove duplicates**
 - Avoids data leakage, class imbalance, and wasted computation
- **Use of Transfer Learning**
 - MobileNetV2 is lightweight and efficient using a pre-trained model for better performance of medical imaging
- **Rotation Invariance Strategies**
 - Add rotation invariance strategies or augmentation tricks for rotated duplicates
- **Wise choice of filename** for better reproducibility, debugging, and interpretation

Thank you!

